# Troponin T hs STAT

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## Kaynaklar

- 1 European patent 394816 and US patent 6376206 by Roche Diagnostics GmbH. Specific antibodies to Troponin T, their production and use in a reagent for the determination of myocardial necrosis.
- 2 Katus HA, Remppis A, Looser S, et al. Enzyme linked immunoassay of cardiac troponin T for the detection of acute myocardial infarction in patients. Mol Cell Cardiol 1989;21(12):1349-1353.
- 3 Katus HA, Scheffold T, Remppis A, et al. Proteins of the troponin complex. Laboratory Medicine 1992;23(5):311-317.
- 4 Alpert JS, Thygesen K, Jaffe A, et al. The universal definition of myocardial infarction: a consensus document: Ischaemic heart disease. Heart 2008;94:1335-1341.
- 5 Hamm CW, Ravkilde J, Gerhardt W, et al. The prognostic value of serum troponin T in unstable angina. N Engl J Med 1992;327(3):146-150.
- 6 Ohmann EM, Armstrong PW, Christenson RH, et al. Cardiac troponin T levels for risk stratification with admission cardiac troponin T levels in acute myocardial ischemia. GUSTO IIA Investigators. N Engl J Med 1996;335:1333-1334.

# Troponin T hs STAT

Troponin T hs (high sensitive) STAT (Short Turn Around Time)

cobas®

- 7 Christenson RH, Duh SH, Newby LK, et al. Cardiac troponin T and cardiac troponin I: relative values in short-term risk stratification of patients with acute coronary syndromes. *Clin Chem* 1998;44(3):494-501.
- 8 Lindahl B, Diderholm E, Lagerqvist B, et al. Mechanisms behind the prognostic value of Troponin T in unstable coronary artery disease: a FRISC II substudy. *J Am Coll Card* 2001;38:979-986.
- 9 Aviles RJ, Askari AT, Lindahl B, et al. Troponin T levels in patients with acute coronary syndromes, with or without renal dysfunction. *N Engl J Med* 2002;346:2047-2052.
- 10 Lindahl B, Venge P, Wallentin L. Troponin T Identifies Patients With Unstable Coronary Artery Disease Who Benefit From Long-Term Antithrombotic Protection. *J Am Coll Cardiol* 1997;29(1):43-48.
- 11 Hamm CW, Heeschen C, Goldmann B, et al. Benefit of abciximab in patients with refractory unstable angina in relation to serum troponin T levels. *N Engl J Med* 1999;340(21):1623-1629.
- 12 Heeschen C, Hamm CW, Goldmann BU, et al. for PRISM study investigators. Troponin concentrations for stratification of patients with acute coronary syndromes in relation to therapeutic efficacy of tirofiban. *Lancet* 1999;354:1757-1762.
- 13 Lindahl B, Diderholm E, Lagerqvist B, et al. Effects on mortality of long-term treatment with l.m.w. heparin in relation to troponin T level and ECG findings - a FRISC 2 substudy. *Eur Heart J* 2000;21(Suppl.):521.
- 14 Newby LK, Ohman EM, Christenson RH, et al. Benefit of Glycoprotein IIb/IIIa Inhibition in Patients With Acute Coronary Syndromes and Troponin T-Positive Status: The PARAGON-B Troponin T Substudy. *Circulation* 2001;103:2891-2896.
- 15 The Task Force for the diagnosis and Treatment of Non-ST-Segment Elevation Acute Coronary Syndromes of the European Society of Cardiology. Guidelines for the diagnosis and treatment of non-ST elevation acute coronary syndromes. *European Heart Journal* 2007;28:1598-1660.
- 16 Thygesen K, Alpert JS, White HD. On behalf of the Joint ESC/ACCF/AHA/WHF Task Force for Redefinition of Myocardial Infarction. *JACC* 2007;50:2173-95.
- 17 Latini R, Masson S, Anand IS, et al. Prognostic Value of Very Low Plasma Concentrations of Troponin T in Patients with Stable Chronic Heart Failure. *Circulation* 2007;116:1242-1249.
- 18 Omland T, De Lemos JA, Christophi C, et al. Distribution and determinants of very low levels of cardiac troponin T in patients with stable coronary artery disease: The PEACE trial. *Eur Heart J* 2008;9(202);1342.
- 19 Omland T, De Lemos JA, Christophi C, et al. Very low cardiac troponin T concentrations and cardiovascular events in patients with stable coronary artery disease: The PEACE trial. *Eur Heart J* 2008;29(202);1644.
- 20 Allen J, Kleiman NS, Nassif D, et al. Prevalence and prognostic significance of preprocedural cardiac troponin elevation among patients with stable coronary artery disease undergoing percutaneous coronary intervention. *Circulation* 2008;118:632-638.
- 21 Sato Y, Kataoka K, Matsumori A, et al. Measuring serum aminoterminal type III procollagen peptide, 7S domain of type IV collagen, and cardiac troponin T in patients with idiopathic dilated cardiomyopathy and secondary cardiomyopathy. *Heart* 1997;78:505-508.
- 22 Dierkes J, Domrose U, Westphal S, et al. Cardiac troponin T predicts mortality in patients with renal disease. *Circulation* 2000;102:1964-1969.
- 23 Ooi DS, Zimmerman D, Graham J, et al. Cardiac troponin T predicts long-term outcomes in hemodialysis patients. *Clin Chem* 2001;47:412-417.
- 24 Apple FS, Wu AHB. Myocardial infarction redefined: Role of cardiac troponin Testing. *Clin Chem* 2001;47:377-379.
- 25 deFilippi C, Wasserman S, Rosanio S, et al. Cardiac Troponin T and C-reactive protein for predicting prognosis, coronary atherosclerosis, and cardiomyopathy in patients undergoing long term hemodialysis. *JAMA* 2003;290:353-359.
- 26 Iliou MC, Fumeron C, Benoit MO, et al. Prognostic value of cardiac Markers in ESRD: Chronic Hemodialysis and New Cardiac markers Evaluation (CHANCE) study. *Am J Kidney Dis* 2003;42:513-523.
- 27 Scott B, Deman A, Peeters P, et al. Cardiac troponin T and malondialdehyde modified plasma lipids in haemodialysis patients. *Nephrol Dial Transplant* 2003;18:737-742.
- 28 Madsen LH, Ladefoged S, Hildebrandt P, et al. Comparison of four different cardiac troponin assays in patients with end-stage renal disease on chronic haemodialysis. *Acute Card Care* 2008;10(3)173-180.
- 29 ver Elst KM, Spapen HD, Nguyen DN, et al. Cardiac troponin I and T are biological markers of left ventricular dysfunction in septic shock. *Clin Chem* 2000;46:650-657.
- 30 Rossing P, Jorsal A, Tarnow L, et al. Plasma hs-Troponin T predicts cardiovascular and all cause mortality as well as deterioration in kidney function in type 1 diabetic patients with nephropathy. Abstract EASD 2008.
- 31 European patent application 890154 by F. Hoffmann-La Roche AG and Roche Diagnostics GmbH. Cardiac Troponin as an indicator of advanced coronary artery disease.
- 32 European patent application 1837659 by F. Hoffmann-La Roche AG and Roche Diagnostics GmbH. Means and methods for the differentiation of acute and chronic myocardial necrosis in symptomatic patients.
- 33 Masson S, Pirelli S, Pulitano G, et al. Circulating Markers of myocyte injury predict first recurrence of atrial fibrillation. *AHA* 2008 November 8-12. New Orleans.
- 34 Lauer B, Niederau C, Kühl U, et al. Cardiac troponin T in patients with clinically suspected myocarditis. *JACC* 1997;30:1354-1359.
- 35 Swaanenburg JC, Klaase JM, DeJongste M, et al. Troponin I, troponin T, CK-MB-activity and CK-MB mass as markers for the detection of myocardial contusion in patients who experienced blunt trauma. *Clin Chim Acta* 1998;272:171-181.
- 36 Giannitsis E, Müller-Bardorff M, Kurowski V, et al. Independent prognostic value of cardiac troponin T in patients with confirmed pulmonary embolism. *Circulation* 2000;102:211-217.
- 37 Herman EH, Lipshultz SE, Rifai N, et al. Use of cardiac troponin T levels as an indicator of doxorubicin-induced cardiotoxicity. *Cancer Res* 1998;58:195-197.
- 38 European patent 1363128 by F Hoffmann-La Roche AG and Roche Diagnostics GmbH. Diagnosis of myocardial infarction and acute coronary syndrome by combination of markers.
- 39 European patent application 1925943 by F. Hoffmann-La Roche AG and Roche Diagnostics GmbH. Means and methods for optimization of diagnostic and therapeutic approaches in chronic artery disease based on the detection of Troponin T and NT-pro BNP.
- 40 Davis GK, Labugger R, Van Eyk JE, et al. Cardiac troponin T is not detected in western blots of diseased renal tissue. *Clin Chem* 2001;47:782-783.
- 41 Ricchiuti V, Voss EM, Ney A, et al. Cardiac Troponin T isoforms expressed in renal diseased skeletal muscle will not cause false positive results by the second generation cTnT assay by Boehringer Mannheim. *Clin Chem* 1998;44(9):1919-1924.
- 42 Hallermayer K, Klenner D, Vogel R. Use of recombinant human cardiac troponin T for standardization of third generation troponin T methods. *Scand J Clin Invest* 1999;59(Suppl 230):128-131.
- 43 World Health Organization. Report of the Joint International Society and Federation of Cardiology/World Health Organization Task Force on Standardization of Clinical Nomenclature. Nomenclature and criteria for diagnosis of ischemic heart disease. *Circulation* 1979;59:607-609.
- 44 Müller-Bardorff M, Hallermayer K, Schröder A, et al. Improved troponin T ELISA specific for cardiac troponin T isoform: assay development and analytical validation. *Clin Chem* 1997;43(3):458-466.
- 45 Kampmann M, Rauscher T, Müller-Bardorff M, et al. Clinical Evaluation of Troponin T and CK-MB Mass on the Elecsys 2010 analyzer. Poster presented at Medlab 97, 12th IFCC European Congress of Clinical Chemistry, August 17-22, 1997, Basel, Switzerland.

# Troponin T hs STAT

Troponin T hs (high sensitive) STAT (Short Turn Around Time)

cobas®

- 46 Mendis S, Thygesen K, Kuulasmaa K, et al. World Health Organization definition of myocardial infarction: 2008-09 revision. *Int J Epidemiol* 2011;40(1):139-146.
- 47 Antmann E, Bassand J-P, Klein W, et al. Myocardial infarction redefined - a consensus document of the Joint European Society of Cardiology/American College of Cardiology committee for the redefinition of myocardial infarction: The Joint European Society of Cardiology/American College of Cardiology Committee. *Am Coll Cardiol* 2000;36:959-969.
- 48 Kobayashii S, Tanaka M, Tamura N, et al. Serum cardiac troponin T in polymyositis/dermatomyositis. *Lancet* 1992;340(8821):726.
- 49 Mair J, Dienstl F, Puschendorf B. Cardiac Troponin T in the Diagnosis of Myocardial Injury. *Critical Reviews in Clin Lab Sci* 1992;29(1):31-57.
- 50 Carrier M, Solymoss BC, Raymond C, et al. Cardiac Troponin T and Creatine Kinase MB Isoenzyme as Biochemical Markers of Ischemia after Heart Preservation and Transplantation. *J Heart Lung Transplant* 1994;13(4):696-700.
- 51 Löfberg M, Tähtelä R, Härkönen M, et al. Myosin Heavy-Chain Fragments and Cardiac Troponins in the Serum in Rhabdomyolysis. Diagnostic Specificity of New Biochemical Markers. *Arch Neurol* 1995;52:1210-1214.
- 52 Anderson JR, Hossein-Nia M, Brown P, et al. Donor Cardiac Troponin T Predicts Subsequent Inotrope Requirements Following Cardiac Transplantation. *Brief Communications in Transplant* 1994;58(9):1056-1057.
- 53 Franz WM, Remppis A, Kandolf R, et al. Serum Troponin T: Diagnostic Marker for Acute Myocarditis. *Letter to the Editor in Clin Chem* 1996;42(2):340-341.
- 54 Mair P, Mair J, Koller J, et al. Cardiac troponin T in the diagnosis of heart contusion. *Lancet* 1991;338:693.
- 55 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem* 1988 Nov;26(11):783-790.

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© 2013, Roche Diagnostics
Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim