



GERMAN-MADE DIAPASONS - TUNING FORKS



ENT TUNING FORKS:

All made in Germany, with age-long craftsmanship experience, in special alloy that ensures maximum sound precision.


GIMA code	ENT TUNING FORKS
31230	Tuning fork 128 Hz* - with weight
31231	Tuning fork 256 Hz* - with weight
31232	Tuning fork 512 Hz
31233	Tuning fork 1024 Hz
31234	Tuning fork 2048 Hz
31235	Tuning fork 4096 Hz
31240	Set of 5 tuning forks in a stylish lined case (contains 128 Hz*, 256 Hz*, 512 Hz, 1024 Hz, 2048 Hz). <i>Tuning forks with asterisk * are supplied with weight</i>

ENT/NEUROLOGY TUNING FORKS

All made in a special alloy and supplied with base.

GIMA code	ENT/NEUROLOGY TUNING FORKS
31245	Tuning fork a1-440 Hz - ENT - 90 g
31246	Tuning fork - LUCAE - 80 g adjustable from c128 to h240 with foot - 120 g
31247	Tuning fork C64, c128 - RYDEL SEIFFER detachable slider either C64 or c128 - 120 g
31249	Hammer neurotone 5 in a lined case

Neurotone hammer: set of 5 components for a complete neurological examination: percussion hammer + tuning fork with weight + pin + brush. All fitted in a stylish case.




ALUMINIUM ALLOY DIAPASONS - TUNING FORKS

- 31225 ALUMINIUM TUNING FORK 128 Hz
- 31226 ALUMINIUM TUNING FORK 128 Hz - with foot

Made from finest aluminium alloy and precisely machine to meet utmost requisites of medical profession.

Resting pad for forehead



NYSTAGMUS SPECTACLES - DR. FRENZEL

• 31204 NYSTAGMUS SPECTACLES - DR. FRENZEL

Easy to use in the search of spontaneous or provoked nystagmus. Fixed, biconvex lens + 18.0 diopters.

Nystagmus spectacles should be connected only with < 4.2V transformer as from safety norms 93/42/EEC, see "GIMA6" at page 166. Made in Germany.

• 31206 BATTERY NYSTAGMUS SPECTACLES

This model allows wireless working and is not bound to power supply. The Nystagmus Spectacle is provided with power through the adjustable battery placed above the spectacle and requires three 1.5 Volt batteries, size C LR14/AM-2 (baby cells). Illumination can be adjusted through a turn of the switch.

• 31130 AC/DC TRANSFORMER - see at page 166

• 31208 NYSTAGMUS BULB - spare

4 Volt 0.230 Ampere, 1 Watt threading Ø 3.5x0.35 mm



TRANSFORMER NOT REQUIRED

STANDARD AND USB VIDEONYSTAGMOSCOPE SYSTEMS

• 31210 VIDEO-NYSTAGMOSCOPE SYSTEM - Spectacles and 1 camera


• 31213 VIDEO-NYSTAGMOSCOPE - USB - mask with Ir camera and USB cable

They facilitate to view and register the Nystagmus through standard vestibular tests (both position and instrumental proof).

These systems allow to view on the monitor or PC/Notebook (31213 only) images of patient eyes during tests and to make Visual Suppression Test (VST) through one red LED light inside the mask

to allow fixation. Innovative technology allows to get quality images. Power 220V (12V DC).

• 31211 ADDITIONAL CAMERA for 31210 To see both eyes on monitor at same time.



752x582 pixels resolution
520 TV lines